

SKID-OK

Ultralight evacuation chair



SKID-OK is an evacuation chair ideal for transfers on descending stairs thanks to the handle on the rear side and slide system with belts. The four wheels also allow the transport on the plane surface.

Specific features

PVC sheet easily cleanable equipped with quick release buckles to adjust the tension
Footrest
Front wheels pivoting with brake
Rear wheels of large diameter to increase the sliding
Mobile headrest sheet to adapt to the patient's height
Reduced storage dimensions

Reference standard

N.A.

Standard equipment

ST00427A	Black pair of belts for chair
----------	-------------------------------

Optional accessories

SK21004E	Storage/carrying bag for SKID-OK
SK14001B	Wall bracket for skid ok

Technical data

Height with extracted handle	1540 ± 10mm
Height with handle closed	1040 ± 10mm
Length	900 ± 10mm
Width	520 ± 10mm
Thickness folded	175 ± 10mm
Weight	10 ± 1 kg
Loading capacity	150 kg
Material	Al, PVC, PE,PU
Front wheels	ø100 mm
Rear wheels	ø 150mm

Class I Med. Dev. according to Dir. 93/42/CEE

SKID OK

SK20003E

CND Classification

Registration number



V0899

349445

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it