

## Skid-E MAX

### Evacuation chair



Skid – E Max is an evacuation chair ideal for transfers on stairs thanks to the integrated slides that allow easy movement on stairs. 4 wheels allow movement on flat surfaces.

#### Specific Features

- PVC Sheet easy to clean
- Slide allow easy movement and control going down and upstairs
- Footrest
- Front wheels pivoting with brake
- Large diameter rear wheels increase sliding action
- Headrest sheet is movable making the device suitable for patient of different heights

#### Accessories

- |          |                               |
|----------|-------------------------------|
| SK21000E | Skid carry bag                |
| SK14000B | WALL BRACKET FOR SKID/PROSKID |

#### Technical data

Height opened	1580 mm
Height folded	1060 mm
Width	530 mm
Length	900 mm
Thickness folded	330 mm
Weight	12,7 kg
Loading capacity	250 kg
Materials	Al, PVC, PE,PU
Front wheels	ø100 mm
Rear wheels	ø 200 mm

Dimensions and weight are subjected to tolerances of ±10mm and ± 1kg

#### Class I Med.Device according to Dir. 93/42/CEE

SKID-E MAX SILVER FRAME  
BLACK SEAT  
SK20014E



CND Classification

V0899

Rev.0 (28/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)