

EVA E5

Evacuation chair



EVA E5 is a chair for transport of patients from the rescue place to the emergency vehicle or to another device intended for ambulance transport.

Specific features

- The foot operated brake applied to the rear wheels, is easy and quick to operate
- Front telescopic handles and foldable and lockable rear handles make it easy to lift
- Upper handle with integrated headrest and belt for head stabilization
- The integrated slide opens easily allowing easy transport on stairs
- The automatic locking system significantly shortens the storage time
- Footrest and armrests to increase patient comfort
- Belts with integrated rewinding system
- The integrated carrying handle allows to transport the device in closed position, helping to reach the rescue place faster
- Seat and backrest in ABS allows quick and easy cleaning
- Seat with slots for additional belts
- Self positioning front wheels
- Optional EN 1789 fastener

Technical data

Length opened w/handles closed	910 ± 10 mm
Length opened w/handles opened	1325 ± 10 mm
Length w/handles and slides opened	1400 ± 10 mm
Width w/rear handles closed	550 ± 10 mm
Height w/upper handle closed	950 ± 10 mm
Height w/upper handle opened	1365 ± 10 mm
Seat height	500 ± 10 mm
Seat and backrest width	400 ± 10 mm
Thickness folded	260 ± 10 mm
Rear wheels	150 mm
Front wheels	100 mm
Materials	Al, Steel, Nylon
Weight	19 ± 0,5 kg
Loading capacity	170 kg

Optional accessories

SK30001E	E-MAX 10G fastener
SK30002E	E-MAX 20G fastener

Class I Med.Dev. according to Dir. 93/42/CEE

EVA E5 EVACUATION CHAIR

SK30009E



CND Classification

V0899

Registration number

852095

Rev.2 (25/05/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it