

## FLASH 07

### Immobilization/extrication Kit



Immobilization kit including SED extrication device, cervical collars, head immobilizer, carry bag and spider straps to be used with a Spinal board.

#### Specific features

- Lightweight and compact extrication device easy to clean
- Head immobilizer with holes to allow ear inspection
- One piece cervical collars with access for tracheotomy
- PVC Carry bag
- Adjustable spider straps to allow easy application to the spinal board

#### Technical data

SED dimensions	900x830 ± 10mm
SED loading capacity	230 kg
Included cervical collars	Sizes Pediatric, Neckless, Regular
Materials	PVC, Nylon, PE
Overall weight	6,1 ± 0,5 kg

Med.dev. CE marked according to dir. 93/42/CE

FLASH 07 - IMMOBILIZATION KIT

WITH SED

SR00103K



Rev.0 (23/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)