

LIFE GUARD

Professional life saving buoy





Rescue buoy to be used as buoyancy aid of victim or rescuer during water rescue.

Specific features		
The tapered hydrodynamic shape, allow easy movements in water		
3 handles to allow easy and firm grip		
The width of the shoulder strap tape, provide more comfort to the rescuer		
Floating and resistant rope		
The color confers great visibility		
The wide shoulder strap doesn't interfere with rescuer's movements and is comfortable		

Technical data		
Dimensions	655x235x120 ± 5 mm	
Rope length	2820 ± 20 mm	
Shoulder strap width	50 mm	
Shoulder strap circumference (not adjustable)	69 cm	
Materials	PE, PP	
Floating body weight	1,2 ± 0,1 kg	
Overall Weight	1,36 ± 0,1 kg	
Maximum buoyancy force	65 N if completely submerged	

Class I med. Dev. According to Dir. 93/42/CEE		
LIFE GUARD ORANGE -		
PROFESS. LIFE-SAVING BUOY	((
SR10000A		
CND Classification	V9099	
Registration number	1163338	