

## LIFE GUARD

### Professional life saving buoy



Rescue buoy to be used as buoyancy aid of victim or rescuer during water rescue.

#### Specific features

The tapered hydrodynamic shape, allow easy movements in water

3 handles to allow easy and firm grip

The width of the shoulder strap tape, provide more comfort to the rescuer

Floating and resistant rope

The color confers great visibility

The wide shoulder strap doesn't interfere with rescuer's movements and is comfortable

#### Technical data

Dimensions	655x235x120 ± 5 mm
Rope length	2820 ± 20 mm
Shoulder strap width	50 mm
Shoulder strap circumference (not adjustable)	69 cm
Materials	PE, PP
Floating body weight	1,2 ± 0,1 kg
Overall Weight	1,36 ± 0,1 kg
Maximum buoyancy force	65 N <sup>if completely submerged</sup>

Class I med. Dev. According to Dir. 93/42/CEE

LIFE GUARD YELLOW -  
PROFESS. LIFE-SAVING BUOY  
SR10001A



CND Classification V9099

Registration number 1163340

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)