

# SHERPA

## Field stretcher carrier



Sherpa are support trolleys for field stretchers which can also be used for short trips within the equipped area. The Plus version is equipped with a system to keep the device in an inclined position. Compatible with series 100/200/300/Slim/Stock stretchers.

### Specific features

- Robust and lightweight structure
- Large wheels diameter to make transport easier
- Placing the stretcher is fast and easy thanks to the 4 quick-locking fasteners
- Compatible with all Spencer spinal boards and pick-up stretchers and all NATO STANAG emergency stretchers
- Low storage space required

### Technical data

Width when open	640 ± 10 mm
Width when closed	300 ± 10 mm
Length when open	995 ± 10 mm
Length when closed	670 ± 10 mm
Height	950 ± 10 mm
Materials	Al, Steel, Nylon
Weight	19 ± 0.5 kg
Load capacity	160 kg

### Accessories

ST00003C	RAIL BAR FOR SHERPA/SHERPA PLUS
ST00004C	SIDE BAG FOR SHERPA/SHERPA PLUS
ST00005C	TANK HOLDER FOR RAIL BAR

Class I MD compliant with UE Reg. 2017/745

SHERPA MULTIFUNCTION CARRIER

ST00001C



CND Classification	V9099
Registration number	413186

Rev.0 (29/09/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)