

## SHERPA PLUS

Field stretcher carrier with handles and Trendelenburg position



Sherpa Plus are support trolleys for field stretchers which can also be used for short trips within the equipped area. The Plus version is equipped with a system to keep the device in an inclined position. Compatible with series 100/200/300/Slim/Stock stretchers.

### Specific features

- Robust and lightweight structure
- Large wheels diameter to make transport easier
- Mechanism to place the stretcher in Trendelenburg anti-shock position
- Placing the stretcher is fast and easy thanks to the 4 quick-locking fasteners
- Compatible with all Spencer spinal boards and pick-up stretchers and all NATO STANAG emergency stretchers
- Low storage space required

### Technical data

Width when open	640 ± 10 mm
Width when closed	300 ± 10 mm
Length when open	995 ± 10 mm
Length when closed	670 ± 10 mm
Height	950 ± 10 mm
Materials	Al, Steel, Nylon
Weight	20 ± 0.5 kg
Load capacity	160 kg

### Accessories

ST00003C	RAIL BAR FOR SHERPA/SHERPA PLUS
ST00004C	SIDE BAG FOR SHERPA/SHERPA PLUS
ST00005C	TANK HOLDER FOR RAIL BAR

Class I MD compliant with UE Reg. 2017/745

SHERPA PLUS FIELD STRETCHER  
CARRIER

**ST00002C**

**CND Classification** V9099

**Registration number** 413187



NATO STOCK N° 6530152053768

Rev.0 (29/09/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)