

PORTA STOCK

Support frame for emergency stretchers



Porta Stock is a support for field stretchers that allows the stretcher to be kept raised from the ground, facilitating small movements and access to the patient by the doctor.

Available in versions with and without wheels, it is compatible with the Stock series stretchers and is also suitable for Spencer stretchers with a maximum width of 575 mm and a minimum length of 1250 mm (not including handles).

Unsupervised parking is not permitted.

Specific features

Robust and lightweight frame
Wheels to allow movement on flat surfaces and brakes on opposite wheels to make it more stable
Reduced thickness when folded, making it ideal for field use
Designed for use with Spencer stretchers Stock Series, is also suitable for use with Spencer emergency stretchers having maximum width of 575 mm and minimum length (handles excluded) of 1250 mm

Technical data

Length when open	1325 ± 10 mm
Width	640 ± 10 mm
Height when open	750 ± 10 mm
Overall dimensions when folded	975 x 640 x 130 ± 10 mm
Wheels	∅ 150 mm
Weight	11 ± 0,5 kg
Capacity	150 kg
Frame material	Steel

Class I MD compliant with UE Reg. 2017/745

PORTA STOCK WITH WHEELS
ST00023C



CND Classification	V0899
Registration number	1241254

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it