


TGS

Field stretcher with fastener



PRIME series stretchers TGS are devices for lifting and carrying people for use in field rescue situations. Parking is permitted preferably when using the parking stops installed in designated locations within the rescue area.

Specific features

- The provided system locks both the front and rear part of the stretcher
- Polyethylene patient board easy to clean
- Adjustable backrest
- Front and rear telescopic handles
- Wheels have large diameter to increase sliding action
- Brakes on outer wheels to stabilize the device when the vehicle is in motion
- Soft mattress to increase patient comfort
- Included belts with plastic hook for restrain the patient
- Large diameter wheels to increase sliding action

Standard equipment

Stretcher fastener

Technical data

| | |
|---------------------------|----------------------------|
| Stretcher length | 1915 mm |
| Stretcher width | 575 mm |
| Possible anchorage points | 2 with Rear stop accessory |
| Materials | Al, PE, Nylon, Steel |

Class I MD compliant with UE Reg. 2017/745

MOD.T008 STRETCHER
ST00027



Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it