


Locking system

For incubator DRÄEGER TI500



Interface for fixation of Dräger TI500 incubators, suitable for all Spencer Incubator trolleys.

Specific features

- Easy locking and unlocking of the incubator with only one movement
- The position on the stretcher can be adjusted if needed

Technical data

Dimensions	580x590x100 ± 10 mm
Materials	Al, Steel, Nylon
Maximum load	100 kg

Note: The interface is not tested according to EN 1789.

LOCKING INTERFACE FOR
DRÄEGER TI500
ST00033

Rev.0 (26/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it