


## SPENCER 250

Lengthwise compact foldable stretcher



Stretcher free of ground supports intended for transfer of patients when they have to be moved on main stretcher or other carrying devices.

### Specific characteristics

Extremely lightweight and durable structure  
 Very low storage space required  
 Lengthwise foldable  
 Easy to clean sheet  
 The high number of handles make transport easier

### Technical data

Length	1890 ± 10 mm
Length folded	950 ± 10 mm
Width	480 ± 5 mm
Thickness opened	30 mm
Thickness folded	60 ± 5mm
Materials	Al,PVC
Weight	5,5 ± 0,3 kg

### Standard accessories

3 Black belts with plastic buckle

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 250 DESERT - FOLDABLE  
 STRETCHER

ST00037A

CND Classification

Registration number



V9099

Rev.0 (11/04/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)