


S-WINGS

Enlargement platform for stretchers



Accessory that allows to increase the width of the patient surface of the Cross, Crossover and Cinco Mas.

Specific features

- Easy to sanitize PVC mattresses
- 3 integrated retractable belts
- Supplied with transport bag
- High resistance stainless steel structure
- Each support is attached to the stretcher by inserting it on supports applied to the stretcher and fixed by means of pins to allow rapid application and removal

Technical data

Materials	Stainless steel/PVC
Dimensions of each support	22x143 cm
Transport bag dimensions	40x32x150 cm

S-WINGS-ENLARG. PLATF. FOR
STRETCHERS CROSS/CINCO
ST00043

Rev.0 (09/12/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it