

STX 12

Bag for emergency stretchers 110



Bag for the transportation and storage of Spencer foldable emergency stretchers 110.

Specific features

PVC fabric easily washable
 Handles covered by a soft plastic material,
 allowing an easy transportation
 Upper zip closure

Technical data

Dimensions	ø180x2030 ± 20 mm
N° handles	2
Material	PVC
Weight	780 ± 20 g

STX 12 TRANSP BAG FOR STRETCHER SPENCER
 110
ST00112A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it