

STX 13

Bag for emergency stretchers 100



Bag for transportation and storage of Spencer foldable emergency stretchers 100.

Specific features

- PVC fabric easily washable
- Handle covered by a soft plastic material, allowing an easy transportation
- Rear belts with holes for wall hanging
- Upper strap closure

Technical data

Dimensions	180x490x1030 ± 20 mm
N° handles	1
Material	PVC
Weight	1000 ± 20 g

STX 13 TRANSP BAG FOR STRETCHER SPENCER
100
ST00113A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it