

STX 14

Bag for emergency stretchers 122



Bag for the transportation and storage of Spencer foldable emergency stretchers 122.

Specific features	Technical data
PVC fabric easily washable	Dimensions 1070x325 ± 10 mm
Handles covered by a soft plastic material, allowing an easy transportation	N° handles 1
Superior strap closing	Material PVC
	Weight 470 ± 0,02 g

STX 14 SACCA TRASP BARELLINI SPENCER 122
ST00114A

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it