

SPENCER 122

Foldable stretcher



Foldable stretcher with ground supports for patient recovery and transfer to a main stretcher or other transport device.

Specific features

Extremely lightweight and durable structure
 Very low storage space required
 Lengthwise and breadthwise foldable
 Easy to clean sheet

Technical data

Dimensions opened	2010x550x140 ± 10 mm
Dimensions folded	1000x400x170 ± 10 mm
Materials	Al, Acciaio,PVC
Loading capacity	170 kg
Weight	6,5 ± 0,3 kg

Standard equipment

ST00598A STX598 RESTRAINT BELTS 2pcs

Reference standards

NA

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 122 – FOLDABLE
 STRETCHER YELLOW
 ST00120A



CND Classification	V9099
Registration number	172358

Rev.0 (02/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it