

SPENCER 122

Foldable stretcher



Spencer 122 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise and breadthwise foldable
- Easy to clean sheet

Technical data

| | |
|-------------------------|----------------|
| Length when open (mm) | 2010 ± 10 mm |
| Length when folded (mm) | 1000 ± 10 mm |
| Width when open (mm) | 550 ± 5 mm |
| Width when folded (mm) | 400 |
| Height when open (mm) | 140 |
| Height when folded (mm) | 170 ± 5 mm |
| Weight (kg) | 6.5 ± 0.3 kg |
| Capacity (kg) | 170 |
| Gripping points (no.) | 4 |
| Support feet (no.) | 4 |
| Materials | Al, Steel, PVC |

Standard equipment

STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

SPENCER 122 - FOUR FOLD
STRETCHER ORANGE
ST00122A



| | |
|---------------------|--------|
| CND Classification | V9099 |
| Registration number | 172360 |

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
spencer.it