


SPENCER 122

Foldable stretcher



Spencer 122 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features
Extremely lightweight and durable structure
Very low storage space required
Lengthwise and breadthwise foldable
Easy to clean sheet

Technical data	
Length when open (mm)	2010 ± 10 mm
Length when folded (mm)	1000 ± 10 mm
Width when open (mm)	550 ± 5 mm
Width when folded (mm)	400
Height when open (mm)	140
Height when folded (mm)	170 ± 5 mm
Weight (kg)	6.5 ± 0.3 kg
Capacity (kg)	170
Gripping points (no.)	4
Support feet (no.)	4
Materials	Al, Steel, PVC

Standard equipment	
STX 598 Restraint belts 2 pcs	

Class I MD compliant with UE Reg. 2017,	/745
SPENCER 122 - FOUR FOLD	
STRETCHER MILITARY	((
ST00124A	
CND Classification	V9099
Registration number	172362