

SPENCER 128

Four foldable stretcher



Foldable stretcher with for patient recovery and transfer to a main stretcher or other transport device.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Widthwise foldable and twice lengthwise foldable
- Transport handles offer a firm grip
- Closure belt to keep the stretcher folded during storage

Standard equipment

- Carry bag
- Couple of black belts

Reference standard

NA

Technical data

Length	2095 ± 5mm
Width	540 ± 5mm
Dimensions folded	580x210x120 ± 10mm
Ground clearance	70 ± 5mm
Belts length	Max 85cm
Materials	Al, Nylon
Loading capacity	159 kg
Weight	6 ± 0,3 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 128 - ORANGE
STRETCHER FOLDABLE IN 4 PARTS
ST00128

CND Classification



V9099

Distributed by Spencer Italia s.r.l.

Rev.0 (14/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it