


SPENCER 205

Compact folding transport stretcher with wheels



Spencer 205 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise foldable
- Foldable supports
- Easy to clean sheet
- The high number of handles and the wheels make transport easier

Technical data

Length when open (mm)	1890 ± 10 mm
Length when folded (mm)	950 ± 10 mm
Width (mm)	480 ± 5 mm
Thickness when folded (mm)	80 ± 5 mm
Height from the ground	180 ± 5 mm
Weight (kg)	6.2 ± 0.2 kg
Capacity (kg)	170
Gripping points (no.)	8
Wheels (no.)	2
Support feet (no.)	2
Materials	Al, PVC

Standard equipment

STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

SPENCER 205 LOW-PROFILE-EASY-FOLD WHEELED STRETCH

ST00205A

CND Classification V9099

Registration number 168604



Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it