

# SPENCER 281

## Emergency stackable stretcher with headrest



Stretcher for recovery of patients and their transfer to a main stretcher or other transport devices.

### Specific features

- Extremely lightweight and durable structure
- Low storage space required thanks to the possibility to stack many stretchers of the same model
- Backrest adjustable in 4 positions in addition to the horizontal one
- Easy to clean PVC sheet
- The dimensions of the tubular frame allow a firm grip during transfer procedures

### Optional accessories

ST00598A STX598 RESTRAINT BELTS

### Technical data

Length	1885 ± 10 mm
Width	480 ± 5 mm
Ground clearance	175 ± 10 mm
Materials	Al, Steel, PVC
Weight	6 ± 0,5 kg
Loading capacity	170 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 281 STACKABLE  
EMERGENCY STRETCHER YELLOW  
W/ HEADREST  
**ST00281A**



CND Classification	V9099
Registration number	94193

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)