

## Spencer 402/3

Transport chair with 3 wheels



Spencer 402/3 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

### Specific features

- Easy to clean PVC sheets
- Rear handles and armrests are obtained from one metal piece allowing to reduce the weight
- Telescopic front handles to facilitate the handling
- The closing system, allows a considerable reduction of the overall dimensions, also in height, making possible the storage in confined spaces
- Strong and light frame
- Footrest and armrests to improve the stability and patient comfort
- Swiveling front wheel

### Standard equipment

ST00427A PAIR OF BELTS FOR CHAIR 420/425

### Optional accessories

ST42703E WALL FIXATION CERTIFIED 10g

### Technical data

Length opened w/closed handles	740 ± 10 mm
Length opened w/opened handles	1040 ± 10 mm
Width	530 ± 5 mm
Height opened	845 ± 10 mm
Height folded	620 ± 5 mm
Thickness folded	180 ± 5 mm
Rear wheels	∅ 125 mm
Front wheels	∅ 100 mm
Materials	Steel/Al/PVC
Weight	9,5 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 402/3 FOLDABLE  
TRANSPORT CHAIR  
ST00402B



CND Classification	V9099
Registration number	93527

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)