

## Spencer 407

### Transport chair with 4 wheels



Spencer 407 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

#### Specific features

- Easy to clean PVC sheets
- Robust and lightweight structure
- Footrest to improve the stability and patient comfort
- Large rear wheels to reduce vibration during patient transfer
- Swiveling front wheels with brake for rotation and axis
- Quick and intuitive closure system
- Safety system against accidental closure

#### Reference standard

NA

#### Standard equipment

ST10423A | 2 belts for chair

#### Optional accessories

ST42705E | 10g wall fastener for chair 407

#### Technical data

Length	810 ± 10 mm
Width	420 ± 10 mm
Height opened	970 ± 10 mm
Height folded	800 ± 10 mm
Thickness folded	250 ± 10 mm
Rear wheels	∅ 200 mm
Front wheels	∅ 100 mm
Materials	Steel/Al/PVC
Weight	10,8 kg ± 0,3
Loading capacity	150 kg

#### Class I Med.Device according to Dir. 93/42/CEE

SPENCER 407 CHROME FOLDING  
EVAC CHAIR W/4 WHEELS



ST00407A

CND Classification | V9099

Registration number | 93552

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)