

STX423


Couple of belts for transport chairs



Adjustable belts to stabilize the patient during transport procedures.

Specific features

Nylon hook to ensure maximum safety
Suitable for use with chairs 406/407/420/425



Technical data ¹

Length L	Max 1300 mm
Belt width	50 mm
Materials	Nylon, PP
One belt weight	130 ± 10 g

¹ The indicated lengths are subject to tolerances of ± 30 mm

Class I Med.Device according to Dir. 93/42/CEE

PAIR OF BELTS 406/407/420/425

CHAIR

ST00423A

CND Classification

Registration number



V9099

98867

Rev.0 (29/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it