

STX427

Couple of belts for transport chairs



Adjustable belts to stabilize the patient during transport procedures.

Specific features


Nylon hook to ensure maximum safety

Suitable for use with chairs 406/407/420/425

Technical data ¹

Length L	Max 1300 mm
Belt width	50 mm
Materials	Nylon, PP
One belt weight	130 ± 10 g

¹ The indicated lengths are subject to tolerances of ± 30 mm



Class I Med.Device according to Dir. 93/42/CEE

PAIR OF BELTS 420/425 CHAIR

ST00427A

CND Classification



V9099

Registration number

98872

Rev.0 (06/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it