

## FIXO KID

### Paediatric restraint system



Fixo Kid is a restraint system for paediatric patients to be used to obtain a level of restraint that allows the patient to be transported safely when used in accordance with the product in use.

#### Specific features

- Compact and lightweight
- Easy to apply
- Easy to clean
- Supplied with carry bag

#### Technical data

Length when extended (mm)	800 ± 10
Width when extended (mm)	485 ± 5
Soft padding thickness	5 mm
Dimensions when placed in bag (mm)	430 x 330 x 80
Weight	2 ± 0.2 kg
Patient retention belts	5
Device fastening straps	3
Minimum load capacity (kg)	4
Maximum load capacity (kg)	20
Main cloth material	PVC
Belts material	Nylon
Hooks material	Stretcher attachment hooks: steel Patient fastening hook: POM

Class I MD compliant with UE Reg. 2017/745

**ST00480C** – FIXO KID - PEDIATRIC RESTRAINT SYSTEM



Rev.0 (09/11/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)