

DNA WALL

Restraint belts



Belt with integrated retractor and female hook which can be used to hold Spencer products.

Specific features

The integrated retractor allows a quick use of the stored product

Quick release buckle locking system

Easy to install, it must be used together with the male system with retractor or fixed

Technical data ¹

Dimensions	85x50x110mm (with buckle placed like shown in the picture)
Belt length	Max 450 mm
Belt width	50 mm
Materials	Nylon, PP, Steel
Weight	160 ± 10 g

¹ Indicated lengths are subject to tolerances of ± 30 mm

DNA WALL F PLASTIC BUCKLE
ST00486F

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it