

## DNA WALL

### Restraint belts



Belt with integrated retractor and male hook which can be used to hold Spencer products.

#### Specific features

The integrated retractor allows a quick use of the stored product

Quick release buckle locking system

Easy to install, it must be used together with the female system with retractor or fixed

#### Technical data <sup>1</sup>

Dimensions	68x87x60mm (with buckle placed like shown in the picture)
Belt length	Max 650 mm
Belt width	50 mm
Materials	Nylon, PP, Steel
Weight	160 ± 10 g

<sup>1</sup> Indicated lengths are subject to tolerances of ± 30 mm

DNA WALL M PLASTIC BUCKLE  
ST00487B

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)