

DNA WALL FIXED

Restraint belt



Belt with female hook which can be used to hold Spencer products.

Specific features

- Easy to install by means of 3 screws (not supplied)
- Quick release buckle locking system
- It must be used together with the male system with retractor or fixed

Technical data ¹

Dimensions	170x60x15 mm (with buckle placed like shown in the picture)
Belt width	50 mm
Materials	Nylon, PP
Weight	38 ± 5 g

¹ Indicated lengths are subject to tolerances of ± 10 mm

DNA WALL F PLASTIC BUCKLE
ST00488F

Rev.0 (04/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it