

## STX 591

### One piece belt



STX 591 restraint belts are aids used to obtain a degree of restraint that allows patients to be transported safely when used in accordance with the product in use.


#### Specific features

Extremely versatile and easy to use thanks to the single metal buckle with quick release system  
 Characterized by a wide adjustment range  
 Metal buckle to ensure maximum safety

#### Technical data <sup>1</sup>

	Min (mm)	Max (mm)
Length L	350	1850
Belt width	50 mm	
Materials	Nylon, chrome-plated steel, plastic	
Weight	250 g	
Belt type	1 piece	
Hook	Metal	

<sup>1</sup> The measures are the maximum reachable. The effective length can vary depending on the device and the knot done to fix the belt to the stretcher



Class I MD compliant with UE Reg. 2017/745

STX 591 - ONE PIECE METAL BELT

ST00591A

CND Classification



V9099

Registration number

96040

Rev.0 (08/01/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)