


STX 592

Two pieces belt



STX 592 restraint belts are aids used to obtain a degree of restraint that allows patients to be transported safely when used in accordance with the product in use.

| Specific features |
|---|
| Metal hook to ensure maximum safety |
| Suitable for self-loading stretchers, basket stretchers and |
| spinal boards |



| Technical data ¹ | | |
|--|-------------------------------------|----------|
| | Min (mm) | Max (mm) |
| Length L | 600 | 1600 |
| Available length after knot has been done ² | 410 | 1330 |
| Belt width | | 50 mm |
| Materials | Nylon, chrome-plated steel, plastic | |
| Belt weight | 250 ± 10 g | |
| Belt type | 2 pieces | |
| Hook | Metal | |

 $^{^{1}}$ The indicated lengths are subject to tolerances of \pm 30 mm

² Approximate value. It varies depending on which device the belt is applied

| Class I MD compliant with UE Reg. 2017/745 | | |
|--|-------|--|
| STX 592 - 2 PC YELLOW BELT | | |
| METAL BUCKLE | ((| |
| ST00592A | | |
| CND Classification | V9099 | |
| Registration number | 96104 | |