

# HM

## Isothermal blanket



Isothermal blanket to be used to contribute to the thermic stabilization of traumatized patient, suffering from heat strokes or hypothermia.

### Specific features

Reusable

The reinforced mesh makes it particularly resistant and increases the weight of the device facilitating its use even in windy conditions

Eylets at the corners allow multiple fixing possibilities

The different characteristics of the two sides of the blanket, allow the use both for hypothermia and heat strokes

### Technical data

Dimensions opened	2000x1400 mm
Package dimensions	370x270x40 mm
Eyelets	øint 17 mm
Weight	560 ± 10 g
Materials	PE, PET, Al

Class I med. Dev. According to Dir. 93/42/CEE

HM - ISOTHERMAL BLANKET

ST00610C



CND Classification	V9099
Registration number	277414

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)