

HT

Isothermal emergency blanket



The isothermal sheet is a medical device used to contribute to the thermal stabilization of patients, both to maintain heat in case of hypothermia and to reflect it in case of heat stroke.

Specific features

Disposable
The different characteristics of the two sides of the sheet allow to use it both in case of hypothermia and heat stroke

Technical data

Dimensions	220 x 160 cm
Weight	62 g
Dimensions folded sheet	9 x 14 cm approx
Materials	Polyester film with vapour deposited aluminium

Class I MD compliant with UE Reg. 2017/745

HT - ISOTHERMAL SHEET

220x160cm

ST00630A

CND Classification



V9099

Registration number

103930

Rev.0 (14/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it