


STX 706

Grey Body Bag



STX 706 is a device intended for recovery and short transfer of corps.

Specific Characteristics

- Heat sealed waterproof PVC fabric
- Handles integrated into the device
- The long zip closure facilitates positioning of corps
- Transparent window for the insertion of identification labels

Technical data

Length	1980 mm
Width	900 mm
N° of handles	6
Material	PVC
Loading capacity	150 kg
Weight	2,6 kg

STX 706 - NYLON GREY BODY BAG
6 HANDLES 89X198CM
ST00706A

Rev.1 (15/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
mortuary.spencer.it