

TRY 763

Pediatric manikin



Paediatric manikin studied to allow an effective training in cardiopulmonary resuscitation (CPR) and external cardiac massage (MCE). The electronic console is equipped with LEDs to indicate if the resuscitation procedures are carried out properly.

Specific Characteristics

Possibility to fill with water to simulate the weight and feel of a real baby

Realistic chest rising

Realistic anatomical areas (sternum, rib cage)

Console with Led for:

Correct finger position

Correct cardiac compression

Correct ventilation

Lever to exclude acoustic alarm

Console compliant with directive EMC 89/336/EEC

Equipment

Carry bag

Clothes

12 spare lungs

Inflation bulb

4 Batteries type AA included

Technical data

Dimensions 82x26x25 cm

Weight 2,3 kg

Console compliant with EN60601-1-2

Console specifications

EN6100-4-2 6 kV Contact

EN6100-4-2 8 kV Air

EN61000-4-4 5kHz EFT/Burst

TRY 763 - PAEDIATRIC MANIKIN

FOR CPR

ST01300C

Distributed by Spencer Italia s.r.l

Rev.1 (24/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it