

TRY 753

Infant CPR manikin



Full body infant manikin studied to allow an effective training in cardiopulmonary resuscitation (CPR) and external cardiac massage (MCE). The electronic console is equipped with LEDs to indicate if the resuscitation procedures are carried out properly.

Specific Characteristics

- Air inflated
- Possibility to fill with water to simulate the weight and feel of a real baby
- Console with Led for:
 - Correct finger position
 - Correct cardiac compression
 - Correct ventilation
- Console compliant with directive EMC 89/336/EEC
- Lever to exclude acoustic alarm

Equipment

- Carry bag
- Clothes
- 12 spare lungs
- Inflation bulb
- 4 Batteries type AA included

Technical data

Dimensions	50x20x12 cm
Wiehgt if air inflated	900 g
Weight if water inflated	2,5 kg
Console compliant with EN60601-1-2	
Console specifications	
EN6100-4-2	6 kV Contact
EN6100-4-2	8 kV Air
EN61000-4-4	5kHz EFT/Burst

**TRY 753 INFANT MANIKIN FOR
CPR WITH ELECTRONIC PAD
AND TRANSPORT BAG
ST01310C**

Distributed by Spencer Italia s.r.l

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it