

TRY 773

CPR mannequin for water rescue



Mannequin designed to provide effective training to CPR procedures in water rescue conditions.

Specific features

The skeleton of the mannequin is made from robust coated steel cables

The body is articulated and allows the use in wider rescue scenarios.

The holes in the limbs, allow water to enter inside increasing the weight of the dummy giving it a more realistic position in water

The chest can be filled with water by opening the valve on the shoulder

If completely filled with water, the manikin will sink to the mouth level

Is possible to apply additional weight on the dummy's chest (not supplied) in order to increase sinking

Supplied with 5 mouth-nose pieces and 5 spare airways

Technical data

Dimensions	165x60x25 ± 20cm
Materials	Vinyl, PE, Steel
Weight empty	About 18 Kg
Weight filled	About 48 kg
Weight needed to sink the mannequin	About 5 kg

TRY 773 - WATER RESCUE MANIKIN

ST01910C

Distributed by Spencer Italia srl

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it