

## FLASH 04 Immobilization Kit



Immobilization kit complete of aluminum and foldable spinal board, cervical collars, disposable head immobilizer and carry bag.

### Specific features

Lightweight and compact spinal board characterized by reduced thickness also when folded for storage  
 Disposable inflatable head immobilizer  
 One piece cervical collars with access for tracheotomy  
 PVC Carry bag with metal rings for wall hanging  
 Smooth easy to clean surfaces

### Technical data

Spinal board dimensions	1880x480x60 ± 10mm
Spinal board dimensions folded	940x480x62 ± 10mm
Spinal board weight	8 ± 0,3 kg
Loading capacity	150 kg
Included cervical collars	Sizes Pediatric, Neckless, Regular

Med.dev. CE marked according to dir. 93/42/CE

FLASH 04 - IMMOBILIZATION KIT  
 WITH PORTA LOC  
 ST02002K



Rev.0 (23/08/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)