

## FLASH 02 Immobilization Kit



Complete immobilization system. By using together spinal board, spider straps and head immobilizer included in this kit, is possible to stabilize and immobilize the patient during rescue operations and ambulance transport.

### Specific Characteristics

Lightweight and compact spinal board
Spider straps with strap closure
Head immobilizer with opening for ear inspection
Soft head immobilizer pillows to increase comfort
Adaptable adult/child head immobilizer base by removing the central padding
Smooth easy to clean surfaces

### Kit components

ST02012B	Rock pin – Orange spinal board with pin
ST02035A	Rock straps – Univ. spider straps for sp.board
SH00200A	Spencer contour – Anatomic head immobilizer

### Accessories

ST02106A	Mark Up – system for permanent personalisation
ST02105A	Fix board - 10 g vertical wall support for sp.board
ST02108A	Fix board - 10 g horiz. wall support for sp.board

### Technical data

Spinal board dimensions	1840x450x50 mm
Spinal board weight	7,8 ± 0,2 kg
Loading capacity	200 kg
Central belt length	1230 mm
Belt weight	460 g
Head block dimensions <sup>1</sup>	235x135x160 mm
Head immob. base dimensions	380x270x20 mm
Head immobilizer weight	1,5 kg
Materials	PE, Nylon, PVC
<sup>1</sup> Size of one headblock	

Med.dev. CE marked according to dir. 93/42/CE

FLASH 02 – SPINAL BOARD ROCK  
PIN+ROCK STRAPS+CONTOUR  
**ST02011K**



Registration number 58063

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)