

PIN STRAPS

Belt system with hooks



Pin Straps full body restraint belts are aids used to immobilise patients on Spencer spine boards.

Specific features

Fastening system with carabiner for quick and easy fastening to spine boards equipped with pins
 Strap closure system to allow the use also with spine boards without pins
 4 transversal belts 3 of which can slide for a better adaptation to device in use

Technical data ¹

Lengths	Min (mm)	Max (mm)
Length of side belts	580 ± 30	970 ± 30
System vertical length	1160 ± 50 mm	
Materials	Polypropylene/Nylon/Metal alloy	
Weight	1500 ± 20 g	
Number of belts	10	
Fastening type	Strap+metal spring catches	

¹ Lengths are referred to the installed device. The real length can vary depending on the device in use with the product.

Class I MD compliant with UE Reg. 2017/745

PIN STRAPS- BELT SYSTEM WITH
HOOKS

ST02015A

CND Classification



V9099

Registration number

99840

Rev.0 (10/02/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it