

## **RSP** Pediatric spider straps system for spine board



RSP full body restraint belts are aids used to immobilise patients on Spencer spine boards.

Specific features
Extremely resistant and robust system
Strap closures to make rescue operations faster
Integrated measuring tape to measure the height of the patient, specifically designed to be used with the paediatric spine board Baby-Go
Sliding transversal bands for greater versatility
Reflective insert in the texture of the belt

Technical data	
Central belt length	750 mm
Minimum length "V" belts	420 mm
Maximum length "V" belts	610 mm
Minimum width	520 mm
Maximum width	610 mm
Integrated measuring tape range	1500 mm
Material	Polypropylene/Nylon
Weight	350 g
Number of belts	10
Fastening type	Strap

Class I MD compliant with UE Reg. 20	)17/745
RSP -IMMOBILIZATION BELTS	
PEDIATRIC BACKBOARD	( (
ST02018A	
CND classification	V9099
Registration number	99844