

## RSP

### Pediatric spider straps system for spine board



Pediatric system for the immobilization of pediatric patients to be used on spinal boards with at least 4 handles on each side.

#### Specific characteristics

- Extremely resistant and robust system
- Strap closures to make rescue operations faster
- Integrated measuring tape to measure the height of the patient, specifically designed to be used with the pediatric spinal board Baby-Go
- Sliding transversal bands for greater versatility
- Reflective insert in the texture of the belt

#### Technical data

Central belt length	750 mm
Minimum length "V" belts	420 mm
Maximum length "V" belts	610 mm
Minimum width	520 mm
Maximum width	610 mm
Integrated measuring tape range	1500 mm
Material	Nylon
Weight	460 g

Class I Med.Device according to Dir. 93/42/CEE

ROCK STRAPS PAEDIATRIC RSP -  
IMMOBILIZATION SYSTEM

ST02018A

CND classification V9099

Registration number 99844



Rev.2 (12/01/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)