

REFLEX STRAPS

Spider straps system
for spine board



Retention system for the immobilization of patients to be used on spinal boards with at least 4 handles on each side.

Specific characteristics

- Extremely resistant and robust system
- Strap closures to make rescue operations faster
- Sliding transversal bands for greater versatility
- Reflective insert in the texture of the belt

Technical data

Number of belts	10
Central belt length	1270 mm
Minimum length "V" belts	270 mm
Maximum length "V" belts	540mm
Minimum width	720 mm
Maximum width	1350 mm
Material	Nylon
Weight	700 ± 20 g

Note: Dimensions subjected to tolerances of ± 20 mm

Class I Med.Device according to Dir. 93/42/CEE

ROCK REFLEX FOR SPINE BOARDS
IN ORANGE COLOR

ST02019A

CND classification V9099

Registration number 99414



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it