


## REFLEX STRAPS

Reflective universal spider straps system for spine boards



Reflex Straps full body restraint belts are aids used to immobilise patients on Spencer spine boards.

### Specific features

Reflective inserts integrated into the fabric facilitate rescue in conditions of poor visibility  
 4 transversal belts, 3 of which are movable to guarantee adaptation to the spine board in use  
 Width adjustment range thanks to the strap closures

### Technical data <sup>1</sup>

Lengths	Min (mm)	Max (mm)
"V" tapes	250	520
Transversal tape 1	750	1350
Transversal tape 2	600	1250
Transversal tape 3	550	1110
Transversal tape 4	480	820
Central belt		1200
Materials	Polypropylene/Nylon	
Weight	600 ± 20 g	
Number of belts	10	
Fastening type	Strap	

<sup>1</sup> Indicated lengths refer to the installed device. The effective length may vary depending on the device on which the product is applied.

Class I MD compliant with UE Reg. 2017/745

REFLEX IMMOBILIZING STRAPS  
 FOR ADULT SPINE BOARDS

ST02020A

CND classification V9099

Registration number 99416



Rev.0 (10/02/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)