


## ROCK STRAPS

Spider straps system for spine board



Rock Straps full body restraint belts are aids used to immobilise patients on Spencer spine boards.

### Specific features

- Extremely resistant and robust system
- Strap closures to make rescue operations faster
- Sliding transversal bands for greater versatility

### Technical data

Central belt length	1230 mm
Minimum length "V" belts	270 mm
Maximum length "V" belts	540 mm
Minimum width	720 mm
Maximum width	1350 mm
Material	Polypropylene/Nylon
Weight	600 ± 20 g
Number of belts	10
Fixation type	Strap

Note: Dimensions subjected to tolerances of ± 20 mm

Class I MD compliant with UE Reg. 2017/745

ROCK-STRAPS – IMMOBILIZATION  
BELTS ADULT BACKBOARD  
ST02035A



CND classification	V9099
Registration number	99422

Rev.0 (10/02/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)