T- STRAPS

Spider straps system for spine board



T-Straps full body restraint belts are aids used to immobilise patients on Spencer spine boards.

Specific features

- Extremely resistant and robust system
- Strap closures to make rescue operations faster
- Fixed transversal belts
- Adjustable longitudinal belt

Technical data	
Belt width	50 mm
Central belt length	From 560 to 1470 mm
Minimum length "V" belts	270 mm
Maximum length "V" belts	540 mm
Minimum width	540 mm
Maximum width	1150 mm
Material	Polypropylene/Nylon
Weight	800 ± 20 g
Number of belts	10
Fastening type	Strap

Note: Dimensions subjected to tolerances of ± 30 mm

Class I MD compliant with UE Reg. 2017/745	
T-STRAPS-IMMOBILISATION BELTS	
ADULT BACKBOARD	
ST02039A	
CND classification	V9099
Registration number	99405

Rev.0 (10/02/2021) UNCHECKED COPY – further revisions will be available on <u>http://support.spencer.it</u> Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

© Copyright - This document is property of SPENCER ITALIA S.r.I and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA S.r.I.