


# B-BAK

## Spinal board



The spinal board Spencer B-Bak is a device for recovery, immobilization and transfer to a main stretcher of patients with suspected spinal injuries.

### Specific features

- The number of handles, their placing and their shape offer an easy grip even in difficult conditions
- The tapered shape allows to place the device inside Spencer basket stretchers
- Smooth and rounded surfaces ensure an effective sanitation
- Robust and light structure to ensure high resistance
- Slides on the lower side to minimize sliding resistance during the application
- The low thickness makes the application and storage easier

### Technical data

Dimensions	1840x405x45 mm
Material	PE
N° handles	14
Loading capacity	180 kg
Weight	6,5 ± 0,3 kg

### Optional accessories

ST02105A	Fix board - 10g vert.fastener for sp. boards
ST02108A	Fix board - 10g horiz. Fastener for sp.boards
ST02035A	Rock straps – Universal belt system
SH00200A	Spencer contour – universal head immobilizer

Class I Med.Device according to Dir. 93/42/CEE

B-BAK SPINAL BOARD ORANGE



ST02040B

CND Classification	V0804
Registration number	101235

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)