

FLASH 03

Immobilization Kit



Complete immobilization system. By using together spinal board, spider straps and head immobilizer included in this kit, is possible to stabilize and immobilize the patient during rescue operations and ambulance transport.

Specific Characteristics

Lightweight and compact spinal board
Spider straps with strap closure
Head immobilizer with opening for ear inspection
Soft head immobilizer pillows to increase comfort
Wide strap surfaces make the head block placing easier
Adjustable frontal and chin straps with soft and high-adhesion paddings to optimize the immobilization
Smooth easy to clean surfaces
Ankle immobilizer to increase patient stability
Integrated Mark-Up for a permanent personalization

Technical data

Spinal board dimensions	1840x405x45 ± 5mm
Spinal board weight	6,5 ± 0,3 kg
Handles	14
Loading capacity	180 kg
Belt weight	800 ± 20g
Head block dimensions ¹	255x85x160 mm
Head immob. base dimensions ²	365x260x20 mm
Head immobilizer weight	1,46 ± 0,07 kg
Ankle imm. dimensions	290x160x80 ± 5mm
Ankle immobilizer weight	0,35 ± 0,01 kg
Materials	PE, Nylon, PVC
¹ Size of one head block	
² Buttonhole excluded	

Kit components

ST02040B	B-Back pin – Yellow spinal board with pin
ST02039A	T-Straps – Univ. spider straps for sp.board
ST00110A	Super blue – Universal head immobilizer
SH00151A	FXA PRO – anat. univers.ankle immobilizer

Med.dev. CE marked according to dir. 93/42/CE

FLASH 03 - B-BAK PIN YELLOW + T-STRAPS+ SUPER BLU + FXA PRO
ST02045B



Accessories

ST02105A	Fix board - 10 g vertical wall support for sp.board
ST02108A	Fix board - 10 g horiz. wall support for sp.board

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it