

B-BAK PIN Spinal board with pins



B-bak Pin spine boards are devices designed to lift and immobilize patients with suspected spinal injuries.

Specific features

- The number of handles, their placing and their shape offer an easy grip even in difficult conditions
- The tapered shape allows to place the device inside Spencer basket stretchers
- Smooth and rounded surfaces ensure an effective sanitation Robust and light structure to ensure high resistance
- Slides on the lower side to minimize sliding resistance during the application
- The low thickness makes the application and storage easier Equipped with pins to allow easy fixing of belts with snap hook

Technical data	
Dimensions	1840 x 405 x 45 mm
Materials	PE
Number of handles	14
Number of pins	8
Load capacity	180 kg
Weight	6 kg

Accessories	
ST02105A	Fix board – 10g vert. fastener for spine boards
ST02108A	Fix board – 10g horiz. fastener for spine boards
ST02035A	Rock Straps – Universal belt system
ST00494B	STX 494 – 3 belts with snap hook
ST02106A	Mark-Up - Permanent personalization system
SH00201A	Spencer Contour – universal head immobilizer

Class I MD compliant with UE Reg. 2017/745		
B-BAK PIN ORANGE - SPINE BOARD		
WITH PIN		
ST02046B		
CND Classification	V0804	

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on http://support.spencer.it Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222 www.spencer.it